



Sani-Cloth® HB

For disinfecting non-porous hard surfaces.

The first germicidal wipe that kills the Hepatitis B virus. Kills over 100 other organisms. Effective against Methicillin Resistant Staphylococcus aureus (MRSA). Kills

Vancomycin Resistant Enterococcus Faccalis (VRE). Effective against Escharichia coli (E coli). Kills HIV-1.

- 2071-P-Q08472 PDI Sani-Cloth® HB, Large, 160/tub, 12 tubs/cs
- 2071-P-W85484 PDI Sani-Cloth® HB,X- Large, 65/tub, 6 tubs/cs
- 2071-P-H21082 PDI Sani-Cloth® HB, Individual, Large, 50/bx, 10 bx/cs
- 2071-U22495 PDI Sani-Cloth® HB, Individual, X- Large, 50/bx, 3 bx/cs

PDI Sani-Cloth® Plus

For disinfecting non-porous hard surfaces. Kills a broad spectrum of organisms faster than "Classic" Sani-Cloth.

- 2071-P-Q89072 PDI Sani-Cloth® Plus, Large, 160/tub, 12 tubs/cs
- 2071-P-Q85084 PDI Sani-Cloth® Plus, X- Large, 65/tub, 6 tubs/cs

Super Sani-Cloth®

Germicidal Hard Surface Disinfectant Wipes. For use in hospitals and other critical care areas where control of the hazards of cross-contamination is of prime importance. Use on surfaces and equipment such as stainless steel, formica, glass, carts, baskets, counters, cabinets and other hard non-porous surfaces. Pre-moistened: each wipe is impregnated with a germicidal solution, no mixing. Pre-measured: for effective, convenient application.

Convenient: unique canister design provides the freedom of one-hand dispensing. EPA approved.

- 2071-P-Q55172 Large, 6"x6.75", 160/tub, 12 tubs/cs
- 2071-P-Q86984 X-Large, 11"x12", 65/tub, 6 tubs/cs
- 2071-P-U87295 Individual Foil Packets, 11.5"x11.75", 50/pk, 3 pk/cs

CaviCide®-Surface Disinfectant/Decontaminant Cleaner

CaviCide® is a convenient, ready-to-use, intermediate-level disinfectant which is effective against TB in 5 minutes and HBV in 3 minutes, viruses (hydrophilic and lipophilic), bacteria (including MRSA and VRE) and fungi. It is safe for all areas of the facility including NICU, operating rooms, isolation rooms, patient care areas and laboratories. CaviCide® is safe for use on non-porous surfaces.



- Broad spectrum of kill
- Non-toxic
- Ready-to-use
- Satisfies OSHA's Bloodborne Pathogen Standards
- Compatible with most medical device material

- 4397-13-1002 CaviCide®, 2 oz with Sprayer, 48 btl/cs
- 4397-13-1008 CaviCide®, 8 oz Spray, 12 btl/cs
- 4397-13-1024 CaviCide®, 24 oz Spray, 12 btl/cs
- 4397-13-1000 CaviCide®, Gallon, 4 btl/cs
- 4397-13-1025 CaviCide®, 2.5 Gallon, 2 btl/cs



CaviWipes™ & CaviWipesXL™ Disinfecting Towelettes

A convenient, ready to use, intermediate-level disinfecting wipe available in two sizes to fit all your needs.

- Towelettes stay fully saturated and will not dry out
- Effective against TB in 5 minutes and HBV in 3 minutes
- Ready to use - towel and CaviCide in one
- Towelettes are safe for all non-porous surfaces
- Sturdy and durable, nonwoven towel will not bunch up during use

- 4397-13-1100 CaviWipes, 160 wipes per can, 12 canisters/case
- 4397-13-1150 CaviWipesXL, 66 wipes per can, 12 canisters/case
- 4397-13-1155 CaviWipesXL, single, 50/box, 6 boxes/case
- 4397-13-1175 Canister Wall Bracket, 12/case



KleenAseptic® Germicidal Surface Cleanser

KleenAseptic® is an aerosol cleanser and disinfectant that dispenses as a foam. It is tuberculocidal, bactericidal and fungicidal. KleenAseptic® is

compatible with non-porous medical and dental equipment and environmental surfaces. Its foaming action allows for better penetration. Its non-toxic synergistic formulation has made KleenAseptic the first choice for disinfection in neonatal units since the 1960's. It is an excellent choice for practitioners who prefer an aerosol disinfectant/cleanser. Satisfies OSHA's

Bloodborne Pathogen Standards

4397-13-1016 KleenAseptic®, 16 oz., 12/cs

Burnishine Germicidal Solution

Germicidal Solution kills dangerous germs. Neutral pH with broad spectrum kill. Effective against TB and HIV (AIDS). EPA

registered. Wipe, spray or soak; disinfects in 10 minutes.

1080-S5N3 Germicidal Solution, Quart, 12/cs
1080-S5N Germicidal Solution, Gallon, 4/cs

MetriMist® Natural Aromatic Deodorizer

MetriMist® eliminates offensive odors by cleaning the air, not masking the smell.

MetriMist® may be sprayed directly in the air, on a surface, or added to laundry water to eliminate odors. It is also effective on irritating odors such as hydrogen sulfide, ammonia, and any other chemical compounds. MetriMist®

is biodegradable, non-toxic and non-flammable.

4397-10-1151 MetriMist®, 1 ounce spray, 100 btl/cs
4397-10-1158 MetriMist®, 8 ounce spray, 12 btl/cs



MetriGuard® General Purpose Cleaner and Disinfectant

MetriGuard® is a ready-to-use and multi-purpose tuberculocidal broad spectrum solution. It is safe and compatible for use on the surfaces of non-critical inanimate objects. MetriGuard® was specially designed as a high-quality housekeeping product. Non-toxic.

4397-10-7506 MetriGuard, 32 ounce spray, 12/cs



MadaCide-FD Disinfectant/Cleaner

A Safe and Effective Alcohol Based Broad Spectrum Formula. Has 6 minute TB, 30 Second HIV and 3 minute Hepatitis Kill Times.

4071-7020 32 oz. Spray Bottle, 12/case
4071-7021 One Gallon Bottle, 4/case

MadaWipe-FDW

Fast Drying Wipes for use in Hospitals, Nursing Homes, Critical Care Areas. One of the most powerful and user friendly cleaning products. Contains Isopropyl Alcohol, 2 - Butoxyethanol, cationic surfactant and water. Safe, easy to use, economically priced.

Biodegradable and non-toxic. Greater stability with longer shelf life. No rinsing necessary. Allow to air dry.

4071-7031 MadaWipe-FDW, 6" x 7", 160 Towelettes/Canister, 12/cs



Wet Wipe Dispenser

0927-WB-116 Top-Pull Wet Wipe Dispenser, Powder Coat 6" x 7" x 6" with 5" Dia. Top Hole, 12/cs



Sklar Ready-to-Use Spray Surface Disinfectant

Sklar Disinfectant is a ready-to-use, highly effective cleaner and disinfectant for use in intensive care units and high-risk patient care areas, surgical suites, infant care units, plasma phoresis facilities and other areas where strict environmental control of cross contamination is critical. Kills TB/Hepatitis B & C in 3 minutes

5868-10-1643 24 oz. Bottles, 6/cs
5868-10-1653 1 Gallon Bottles, 4/cs



TimeMist® TLC Tank Flush

- Converts tank toilets to "Touch Free"
- Reduces cross-contamination of germs
- Upgrades washroom image
- Invisible to the patron, vandal resistant
- Manual flush capability
- Automatic flushing provides a cleaner restroom
- Comes complete with all installation tools
- 5 minute installation

6930-35-3701TM TimeMist® TLC Tank Flush – Touchless control tank flush, 6/cs

TimeMist® Commode Cleaning System

The TimeMist Commode Cleaning System utilizes an innovative refill cartridge that cleans and deodorizes both the inside of the tank and the entire commode.



Automatically eliminates odors at the source and keeps the fixture clean. Every month, simply insert the refill cartridge into the commode bracket and hang in the tank. At the end of the month, simply throw away and insert a new one to keep the commode clean and fresh.

6930-35-3265TM TimeMist® Commode Cleaning System, 12/cs

Bodily Fluid Clean Up Pack

One time use pack.
Contains necessary

products to provide personal protection at a bodily fluid spill scene. Pack is small enough to include in a first aid cabinet or store in a convenient location.

9510-214-P Bodily Fluid Clean Up Pack, 16 pc., 1 ea.



EZ-Cleans™ Mercury Spill Kit

In the event of a mercury spill, be prepared to safely contain, clean up, and dispose of the vaporous liquid metal. For use in hospitals, laboratories, clinics, schools or anywhere there is a risk of exposure to a mercury spill. The kit contains Green Z™ to surround the liquid metal, immobilizing and controlling the spill.

Kit contains:

- One 4 oz. Activator Bottle
- Two Pick up Scoops with Scraper
- One Zip Lock Disposable Bag
- One Disposable Bag with Closure
- Two Nitrile Gloves
- One Safety Shield
- 2 oz. Green Z™
- One Towel

6709-P00-19105 EZ-Cleans™ Mercury Spill Kit, 12/cs

The Products Listed In This Catalogue Represent Only A Portion Of Those Available. Please Call Our Office Or See Your Sales Representative.

Ozium® Metered Air Sanitizers

The active ingredients in Ozium® actually cleanse the air and effectively eliminate bacterial and viral pathogens with glycolized action. Lighter-than-air particles are dispersed into the air where they linger for up to 4 hours, encapsulating airborne bacteria and odors. Ozium® effectively sanitizes, deodorizes, and leaves a fresh fragrance in the air. Ozium 3000 metered air sanitizers are for use in the Waterbury TimeMist® Classic and TimeMist® Plus dispensers. Fight airborne bacteria 24/7 for a full 30 days with each can of Ozium® 3000. Available in original, vanilla, fresh citrus, and country fresh.

6930-53-031CWD Ozium® 3000 Original, 12/cs
 6930-53-131CWD Ozium® 3000 Vanilla, 12/cs
 6930-53-231CWD Ozium® 3000 Fresh Citrus, 12/cs
 6930-53-331CWD Ozium® 3000 Country Fresh, 12/cs



TimeMist® Classic Dispenser

Simply reliable. The TimeMist® Classic is a simple and reliable automatic metered dispenser. It provides maximum odor control at a minimum cost. For use with all Ozium 3000 metered refills.

Features:

- Basic on/off control
- 15 minute spray setting
- 30 day refill life
- 5 year guarantee

6930-32-1141TM TimeMist® Classic Dispenser, 12/cs



TimeMist® Plus Dispenser

Sophistication and intelligence in a programmable dispenser. The TimeMist® Plus utilizes cutting-edge electronics to offer simple and complete programmability. It allows you to control fragrance intensity using spray settings (LT, MED, HVY) that fit your odor control needs. Saves time, reduces labor and lower costs. Works with all TimeMist® metered fragrances and Ozium® metered 3000 refills.

Programming Features:

- Time of day
- Start/stop times
- Light, Med. Or Heavy spray settings

Other Features:

- LED service reminder
- Audio service reminder
- LED battery replacement reminder
- 5 year warranty

6930-32-0141TMO4B TimeMist® Plus Dispenser, 12/cs



Ozium® 1500 3.5oz. Hand-held Air Sanitizer

The active ingredients in Ozium® Hand-held 3.5oz sprays cleanse the air with glycolized action. Lighter-than-air particles are dispersed into the air where they linger for up to 4 hours, encapsulating airborne bacteria and odors. Ozium® effectively sanitizes, deodorizes, and leaves a fresh fragrance in the air. Ozium® has been proven effective in biological laboratory testing for hospitals, sick rooms, nursing homes and clinics. Available in original, vanilla, fresh citrus, and country fresh.

6930-53-015CWD Ozium® 1500 original, 36/cs
 6930-53-115CWD Ozium® 1500 vanilla, 36/cs
 6930-53-215CWD Ozium® 1500 fresh citrus, 36/cs
 6930-53-315CWD Ozium® 1500 country fresh, 36/cs



Ozium® 500m 0.8oz. Hand-Held Air Sanitizer

The active ingredients in Ozium® Hand-held 0.8oz. spray cleanse the air with glycolized action. Lighter-than-air particles are dispersed into the air where they linger for up to 4 hours, encapsulating airborne bacteria and odors. Ozium® effectively sanitizes, deodorizes, and leaves a fresh fragrance in the air. Ozium® has been proven effective in biological laboratory testing for hospitals, sick rooms, nursing homes, and clinics.

6930-53-011CWD Ozium® 500m 0.8oz. hand-held spray, 144/cs



Ozium® 7000 14.5oz. Hand-Held Air Sanitizer

The active ingredients in Ozium® Hand-held 14.5oz. spray cleanse the air with glycolized action. Lighter-than-air particles are dispersed into the air where they linger for up to 4 hours, encapsulating airborne bacteria and odors. Ozium® effectively sanitizes, deodorizes, and leaves a fresh fragrance in the air. Ozium® has been proven effective in biological laboratory testing for hospitals, sick rooms, nursing homes, and clinics.

6930-53-017CWD Ozium® 7000 14.5oz. hand-held spray, 12/cs

