

## Relief Pak™ Moist Heat Steam Packs

Simple to use; place pack in heating unit or in a pan of hot water for home treatment. Remove and wrap with a terry cover or regular towel, and apply to treatment area! Store packs in a plastic bag, place in your freezer.

|              |                          |
|--------------|--------------------------|
| 2347-11-1310 | Standard, 10" x 12"      |
| 2347-11-1311 | Neck, 9" x 24"           |
| 2347-11-1205 | Circular, 10" diameter   |
| 2347-11-1312 | Oversize, 15" x 24"      |
| 2347-11-1313 | Halfsize, 5" x 12"       |
| 2347-11-1314 | Spine small, 10" x 18"   |
| 2347-11-1315 | Spine large, 10" x 24"   |
| 2347-11-1316 | Knee-shoulder, 10" x 20" |
| 2347-11-1317 | Cervi-back, 22" x 24"    |
| 2347-11-1318 | Myofacial, 6" x 25"      |
| 2347-11-1319 | OB-GYN (pair), 3" x 11"  |
| 2347-11-1320 | Eye (pair), 3" x 3"      |

## Hydrocollator® Heating Units

Made of stainless steel. Fully automatic, thermostatically controlled. Units maintain packs in water at proper temperature for maximum therapeutic benefits. No special plumbing connections required. 120 volts, 60 Hz. Each unit supplied with indicated number of standard size steam packs below.



|           |   |
|-----------|---|
| 2347-2502 | Mobile Model M-4, 24 Pk Cap, 20" W x 32" L x 33" H      |
| 2347-2402 | Mobile Model M-2, 12 Pk Cap, 15-1/2"W x 23-1/2"L x 31"H |
| 2347-2302 | Mobile Model SS-2, 8 Pk Cap, 15-1/2"W x 21"L x 31"H     |
| 2347-2102 | Stationary Model E-1, 4 Pk Cap, 8" D x 12-1/2"W x 15" H |

## Hydrocollator® Moist Steam Heat Packs and Covers

### Packs

|              |                     |
|--------------|---------------------|
| 2347-00-1006 | Standard, 10" x 12" |
| 2347-00-1002 | Neck, 24"           |
| 2347-00-1004 | Oversize, 15" x 24" |
| 2347-00-1014 | Half Size, 5" x 12" |

### Terry-Filled Covers

|              |                     |
|--------------|---------------------|
| 2347-00-1118 | Standard, 20" x 24" |
| 2347-00-1104 | Neck, 9" x 24"      |
| 2347-00-1126 | Oversize, 24" x 30" |

## COLDSTAR Versatile Cryotherapy Packs



Reusable, Microwavable, Stays Flexible, May be used Hot or Cold

|           |   |
|-----------|---|
| 9811-0801 | Standard Size, 6" x 9", 24/case           |
| 9811-0702 | Junior Size, 4" x 7", 48/case             |
| 9811-0703 | Throat Perineal Size, 4.5" x 10", 24/case |
| 9811-0706 | Ice Bandage, 7.5" x 15", 6/case           |
| 9811-0807 | Ice Bandage, 4.5" x 15", 12/case          |

Note: All sizes are approximate.

## COLDSTAR Reusable Hard Freeze Gel Pack

Used in transport of all perishable items. Ideal for maintaining proper temperature during transport. Special ventilated boxes for quick case freezing.



|             |  |
|-------------|--|
| 9811-170872 | 8oz Transport Gel, 4" x 7", 72/case    |
| 9811-171636 | 16oz Transport Gel, 6" x 8", 36/case   |
| 9811-172424 | 24oz Transport Gel, 6" x 8.5", 24/case |

Note: All sizes are approximate.

## dj Ortho/ProCare Head Halters

|               |   |
|---------------|---|
| 5247-79-83060 | Universal Head Halter, Repro-type head halter. Contoured, padded with adjustable chin piece.                        |
| 5247-79-83050 | Orthopedics Deluxe Head Halter, Durable non-stretch design with adjustable straps. Fully padded, one size fits all. |


## dj Ortho/ProCare Bed Traction Kit

Includes all components for in-bed cervical traction: Plastic 20 pound Water Weight Bag (90132), 12" Spreader Bar with 8' Cord (90134), Padded Head Halter (83041).

|               |                  |
|---------------|------------------|
| 5247-79-90140 | Bed Traction Kit |
|---------------|------------------|

## dj Ortho/ProCare Overdoor Cervical Traction Kit

Includes all items necessary to administer cervical traction at home or in physician's office. Universal kit includes: Plastic 20 pound Water Weight Bag (90132), 12" Spreader Bar with 8' Cord (90134), and Padded Head Halter (83040).



|               |                                |
|---------------|--------------------------------|
| 5247-79-90130 | Overdoor Cervical Traction Kit |
|---------------|--------------------------------|

## dj Ortho/ProCare Full Foam Bucks Traction

Sewn-in medial and lateral stay, 6" spreader bar, loop/lock closure.

- 5247-79-83062 Bucks Traction, XS
- 5247-79-83063 Bucks Traction, S
- 5247-79-83065 Bucks Traction, M
- 5247-79-83067 Bucks Traction, L
- 5427-79-83068 Bucks Traction, XL



## Overdoor Traction Sets

- Provides easy to use cervical traction therapy
- One size foam padded head halter fits most
- Each set includes everything needed for installation: pulleys, traction rope, notched adjustable spreader bar, head halter, S-hook and vinyl weight bag with graduated scale up to 20 lbs.



- 4050-534-2014-0000 Deluxe Overdoor Traction Set with 20 lb. Vinyl Bag and Hook and Loop Closure Head Halter
- 4050-534-2018-0000 20 lb. Universal Vinyl Water Weight Bag
- 4050-534-2020-1900 Hook and Loop Style Closure Head Halter (one size fits most)

## Chester LiquaSonic Ultrasound Gel

LiquaSonic Ultrasound Transmission Gel has a thicker, no-slip formula that remains viscous even when applied to salty skin. Salt free, formaldehyde-free formula eliminates the air between the transducer and skin, offering exceptional coupling efficiency for more effective treatment/diagnosis.

LiquaSonic Gel has been formulated to be acoustically correct for the broad range of frequencies used in medical ultrasound applications. Formula is water soluble and non-greasy preventing damage to expensive transducers and clothing.

- 1365-001205 8 fl. oz. Dispensing Cap, 12/cs
- 1365-001222 5-Liter Collapsible Container, also includes 8 oz. Dispensing Bottle and Cap, 4/cs



## Chester LiquaSonic Ultrasound Lotion

Similar in use to the ultrasound gel, the lotion is often preferred by medical professionals. The aqueous formulation is pleasing to the patient and is acoustically correct for the large spectrum of frequencies used in medical ultrasound applications. The lotion is hypoallergenic, bacteriostatic, salt-free, formaldehyde-free and will not harm clothing, equipment or skin.

- 1365-001268 8 fl. oz. Dispensing Cap, 12/cs
- 1365-001267 128 fl. oz., 4/cs



## Aquasonic® 100 Ultrasound Gel

Aquasonic® 100 is the most widely used gel for diagnostic and therapeutic medical ultrasound. Recommended for all procedures where a viscous gel is required. Hypoallergenic, bacteriostatic and non-sensitizing. Will not stain clothing, damage transducers or irritate skin.

- 4926-01-02 60g Doppler® size tube, 12/bx
- 4926-01-08 0.25 liter dispenser, 12/bx
- 4926-01-34 1 liter with dispenser cap, 6/bx
- 4926-01-50 5 liter SONICPAC® with dispenser, 1/bx, 4bx/cs
- 4926-01-01 20g Overwrapped, sterilized foil pouches, 48/bx

## Scan® Ultrasound Gel

Our exclusive dual-viscosity formula allows SCAN® to be applied as a gel and scan as a fluid. For medical ultrasound procedures where a less viscous gel is desired. Hypoallergenic and bacteriostatic. Will not stain clothing, damage transducers or irritate skin.

- 4926-11-08 8 oz. Dispenser size bottle, 12/bx
- 4926-11-28 1 U.S. gallon with dispenser, 4/cs
- 4926-11-28S SCANPAC® contains: 4 SCAN® gallons, 2 dispenser bottles and 1 dispenser pump